Maintrac[®]

Células Epiteliais Tumerais Circulantes Vivas
(CETCs / CTCs)

[maintrac[®]](http://www.maintrac.com)

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Monitoring

J Clin Oncol 2008, 26 (8):
1208-1215

Monitoring the Response of Circulating Epithelial Tumor Cells (CETC) to Adjuvant Chemotherapy in Breast Cancer Allows Detection of Patients at Risk of Early Relapse

VOLUME 26 · NUMBER 8 · MARCH 10 2008

JOURNAL OF CLINICAL ONCOLOGY

ORIGINAL REPORT

Monitoring the Response of Circulating Epithelial Tumor Cells to Adjuvant Chemotherapy in Breast Cancer Allows Detection of Patients at Risk of Early Relapse

Katharina Pachmann, Oumar Camara, Andreas Kavallaris, Sabine Krause, Nele Malarski, Mieczyslaw Gajda, Tarsten Krall, Cornelia Jo'ne, Ulrike Hammer, Annelore Altendorf-Hofmann, Carola Rabenstein, Ulrich Pachmann, Ingo Runnebaum, and Klaus Ho'fken

A B S T R A C T

From the Clinic for Internal Medicine II, Institution for Pathology, and Women's Hospital, Friedrich Schiller University; Tumorzentrum, Jena; and Transfusionsmedizinisches Zentrum, Bayreuth, Germany.

Submitted July 25, 2007; accepted November 2, 2007.

Authors' disclosures of potential conflicts of interest and author contributions are found at the end of this article.

Corresponding author: Katharina Pachmann, MD, PhD, Department of Experimental Hematology and Oncology, Clinic for Internal Medicine II, Friedrich Schiller Universität Jena, Erlanger Allee 101, D-07747, Jena, Germany; e-mail: katharina.pachmann@med.uni-jena.de.

© 2008 by American Society of Clinical Oncology

0732-183X/08/2608-1208/\$20.00

DOI: 10.1200/JCO.2007.13.6523

Purpose

To demonstrate that it is possible to monitor the response to adjuvant therapy by repeated analysis of circulating epithelial tumor cells (CETCs) and to detect patients early who are at risk of relapse.

Patients and Methods

In 91 nonmetastatic primary breast cancer patients, CETCs were quantified using laser scanning cytometry of anti-epithelial cell adhesion molecule-stained epithelial cells from whole unseparated blood before and during adjuvant chemotherapy.

Results

Numbers of CETCs were analyzed before therapy, before each new cycle, and at the end of chemotherapy. The following three typical patterns of response were observed: (1) decrease in cell numbers (> 10-fold); (2) marginal changes in cell numbers (< 10-fold); and (3) an (sometimes saw-toothed) increase or an initial decrease with subsequent reincrease (> 10-fold) in numbers of CETCs. Twenty relapses (22%) were observed within the accrual time of 40 months, including one of 28 patients from response group 1, five of 30 patients from response group 2, and 14 of 33 patients from response group 3. The difference in relapse-free survival was highly significant for CETC (hazard ratio = 4.407; 95% CI, 1.739 to 9.418; $P < .001$) between patients with decreasing cell numbers and those with marginal changes and between patients with marginal changes and those with an increase of more than 10-fold (linear Cox regression model).

Conclusion

These results show that peripherally circulating tumor cells are influenced by systemic chemotherapy and that an increase (even after initial response to therapy) of 10-fold or more at the end of therapy is a strong predictor of relapse and a surrogate marker for the aggressiveness of the tumor cells.

J Clin Oncol 26:1208-1215. © 2008 by American Society of Clinical Oncology

INTRODUCTION

Solid tumors can seed tumor cells into the peripheral blood, which may, even after complete resection of the tumor, eventually grow into frequent cause of death in women in the developed world. Although early detection, precise surgery with wide margins, and adjuvant therapy have improved results,¹ relapse is not infrequent. In premenopausal women, a first narrow peak with metastatic disease with higher cell numbers in blood.² In metastatic disease, the clinical course, and a second peak occurs at 28 to 30 months. Postmenopausal patients display a peak at approximately 18 to 20 months.² After diagnosis of metastatic disease, the outcome is fatal. To date, there is no tool to monitor the effect of isolated tumor cells in bone marrow apart from statistical analysis, indicating that a portion of the circulating epithelial tumor cells (CETCs) may be biologically irrelevant and tumor cells may differ in

1208

Information downloaded from jco.ascopubs.org and provided by F Delbanco on March 9, 2008 from 141.35.226.179. Copyright © 2008 by the American Society of Clinical Oncology. All rights reserved.

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